Box 4-11 provides a detailed description of these devices.

## •Nursing Alert

The belief that core temperature can be estimated by adding 1° C to the temperature taken in the axilla is incorrect. Do not add a degree to the finding obtained by taking a temperature by the axillary route.

### **Box 4-11**

# Types of Thermometers Used to Measure Temperature in Infants and Children

#### **Electronic Thermometer**

Temperature is sensed with an electronic component called thermistor mounted at the tip of a plastic and stainless steel probe, which is connected to an electronic recorder. A disposable plastic cover is used for infection control.

Temperature measurement appears on digital display within 60 seconds.

The probe can be placed in the mouth, axilla, or rectum.

## **Infrared Thermometer**

Thermal radiation is measured from the axilla, ear canal, or tympanic membrane.

Temperature measurement appears on the digital display in approximately 1 second.

Three types are available for ear-based use: Tympanic, ear canal, and arterial heat balance via the ear canal (AHBE).

Often these devices are all inappropriately referred to as *tympanic thermometers*.

Temperatures measured in this way reflect arterial (bloodstream)